

Light Bond™

Light Cure Orthodontic Adhesive



Contents of Kit

- 20 g Paste • [2] 3 cc Sealant Resin •
- 23 g Liquid Etching Agent •
- Mixing Pad • Spatulas •
- Brushes (Sponge pellets may be substituted for brushes) •

Product Features

- **Available With or Without Fluoride** - Light Bond™ Paste and Sealant are formulated with a patented monomer that contains fluoride. It is available without fluoride in Push Syringe only.
- **Unsurpassed Bond Strength** - The Light Bond™ Paste is a highly filled material that delivers maximum bond strength in all situations with all types of appliances
- **Smooth Tacky Paste Viscosity** - The viscosity of Light Bond™ Paste is very tacky which prevents bracket flotation. If a thinner paste viscosity is desired for bonding lingual retainers and acrylic appliances, Light Bond™ Thin or Medium Paste is available. Specify when ordering.
- **Total Paste Polymerization Under All Bracket Bases** - Light Bond™ is comprised of a catalyst and filler that reflects the light to insure complete polymerization under any appliance: metal, ceramic, sapphire or plastic. Once the material has been properly cured, an arch wire can be placed immediately.
- **Ease of Removal** - Light Bond™ is a small particle glass filled material designed specifically for orthodontic use that allows maximum penetration of filler into the finest screen mesh base and removal without enamel scarring.
- **Suggested Uses** - Bracket adhesive

Light Bond™ Kits

Light Bond™ kit with Paste in Push Syringes with fluoride.....#LBEPF

Light Bond™ kit with Paste in Screw Syringes with fluoride.....#LBESF

Light Bond™ kit with Paste in Tips with fluoride#LBETF

NOTE: Light Bond™ Filled Sealant may be substituted for standard Light Bond™ Sealant. Non-fluoride Release Formula available upon request. Light Bond™ Paste is available in Regular, Medium and Thin viscosities. Light Bond™ Paste can be used with any cured sealant. Light Bond™ Kit is available with Assure® Plus All Surface Bonding Resin, Assure® Universal Bonding Resin (hydrophilic) or PRO SEAL® or L.E.D. PRO SEAL® Enamel Protective Sealant as replacements for Light Bond™ Sealants at an additional charge. Blue shade available upon request.

Individual Light Bond™ Components

5 g Light Bond™ Paste in Push Syringe with fluoride#LBPPF

5 g Light Bond™ Paste in Screw Syringe with fluoride.....#LBPSF

5 g Light Bond™ Paste in Tips with fluoride#LBPTF

3 cc Light Bond™ Sealant with fluoride#LBS3F

7 cc Light Bond™ Sealant with fluoride#LBS7F