


INTRODUCING
A REMARKABLE BREAKTHROUGH
IN CLASS III CORRECTION



INTRODUCING THE ALL-NEW
CARRIERE® MOTION™ CLASS III APPLIANCE

The *Carriere Motion* Class III Appliance provides a new, remarkably easy-to-use and patient-friendly solution for Class III treatment. This discreet, comfortable appliance is direct bonded in just minutes, and is as easily tolerated as elastics alone! Imagine an appliance that gives you and your Class III patients an option without surgery or cumbersome, uncomfortable, and unsightly extra-oral devices. If you've ever struggled while tackling Class III cases, take a look at the all-new *Motion* Class III Appliance today!

Learn more about the *Motion* Class III Appliance at
888.851.0533 or **HenryScheinOrtho.com**.

Pre-Treatment



In Treatment



Post-Treatment



HENRY SCHEIN®
ORTHODONTICS

The New Movement In Orthodontics™

CARRIERE MOTION CLASS III APPLIANCE

The *Motion* Class III Appliance is a minimally invasive appliance intended to treat Class III malocclusions without extractions, painful orthognathic surgery, or traditionally uncomfortable and obtrusive face-masks.* The *Motion* Class III Appliance is as easy to place as a bracket or buccal tube, and is designed to provide predictable results and create a harmonious, balanced profile, all while preserving the patient's natural features. The patient-friendly appliance is used at the beginning of treatment, prior to braces or aligners, and delivers rapid Anterior-Posterior correction to help reduce overall treatment times by several months. Intra-oral elastics connect the appliance with maxillary anchorage to activate the mandibular posterior segment and to move it bodily as a block for an ideal Class I malocclusion. Once the ideal platform is achieved, brackets or aligners are then placed to finish the case.

Designed to provide predictable results and create a harmonious, balanced profile.



KITS		ITEM NUMBER	
<i>Motion</i> Class III Standard Kit (1 set of each 23, 25, 27 mm)		424-407C	
DESCRIPTION/SIZE	COLOR-CODE	ITEM NUMBERS 1 pair	ITEM NUMBERS SINGLE (UNIVERSAL L/R)
<i>Motion</i> Class III 16 mm	White	424-416C	424-416
<i>Motion</i> Class III 18 mm	Blue	424-418C	424-418
<i>Motion</i> Class III 20 mm	Pink	424-420C	424-420
<i>Motion</i> Class III 23 mm	Yellow	424-423C	424-423
<i>Motion</i> Class III 25 mm	Red	424-425C	424-425
<i>Motion</i> Class III 27 mm	Green	424-427C	424-427
ACCESSORIES			
Typodont <i>Motion</i> Class III with Clear Aligner	-	-	631-033DNE
Typodont <i>Motion</i> Class III Maloccluded with Clear Aligner	-	-	631-034DNE
<i>Motion</i> Class III 5X Demonstration model	-	-	631-433X
<i>Motion</i> Elastics - Force 1, 1/4", 6 oz (50 packs of 100 elastics)	-	-	424-9F1
<i>Motion</i> Elastics - Force 2, 3/16", 8 oz (50 packs of 100 elastics)	-	-	424-9F2



For more information on our products and educational offerings, please contact us:
OrthoOrganizers.com
Outside the U.S.: +(1) 760 448 8600 • **Canada:** CERUM 800.661.9567
To email an order outside the U.S.: IntlOrders@OrthoOrganizers.com
 Ortho Organizers, Inc. 1822 Aston Ave. Carlsbad, CA 92008-7306 USA

