


Ultra Band-Lok®

Single Paste Compomer Band Cement
Contains Fluoride



Also available in Luer Loc Syringe



Contents of Ultra Band-Lok® Kit

- 30 g Paste in Push Syringes •

Contents of Luer Loc Ultra Band-Lok® Kit

- 30 g Paste in Luer Loc Push Syringes •
- (60) 15 gauge Luer Loc Dispensing Tips •

Product Features

- **Single Paste System** - No wasteful mixing
- **Light Activated** - No working time constraints
- **No Mixing**
- **All Purpose Cement** - Ultra Band-Lok® can also be used to bond large acrylic appliances and as an occlusal build-up material. When bonding the appliances you must first etch the enamel surface.
- **Bonds to Metal** - Ultra Band-Lok® will bond chemically to stainless steel bands. Removal is all in one piece with cement in the appliance

Available in Natural or Blue Shade

Ultra Band-Lok® Kits

- Ultra Band-Lok® kit in Natural shade with fluoride.....#UBL
 Ultra Band-Lok® kit in Hi Vis Blue shade with fluoride.....#UBLB
 Ultra Band-Lok® Luer Loc w/Dispensing Tips with fluoride.....#LLUBL
 Ultra Band-Lok® Luer Loc Hi Vis Blue shade w/Dispensing Tips with fluoride.....#LLUBLB

Individual Ultra Band-Lok® Components

- 5 g Ultra Band-Lok® Syringe with fluoride.....#UBLPP
 5 g Ultra Band-Lok® Hi Vis Blue shade Syringe with fluoride..#UBLBPP
 5 g Ultra Band-Lok® Luer Loc Syringe w/Dispensing Tips with fluoride..#LLUBLP
 5 g Ultra Band-Lok® Luer Loc Hi Vis Blue shade Syringe w/Dispensing Tips with fluoride.....#LLUBLBP
 (60) 15 gauge Luer Loc Dispensing Tips.....#LLTIPS